

**FACE MASK - 3 PLY SURGICAL FACE MASK - TYPE IIR. PRICE PR. MASK.****\$0.50 MOQ 100.000 masks**

Disposable antibacterial surgical 3 Ply PP Face Mask with earloop for use in factory, house, clinic, hospital, pharmacy, restaurant, food processing, beauty salon, electronics. Available in blue, white, green and pink.

SKU: CHM1**Categories:** [Bioide](#), [FaceProtect](#), [Masks](#), [PersonProtect](#)**PRODUCT DESCRIPTION**

Disposable 3 Ply PP Face Mask With Earloop

- 1). Nose Bar Adaptable
- 2). BFE: 95-99%
- 3). Color:Blue/White/Green/Pink
- 4). Applications:Used in factory, house, clinic, hospital, pharmacy, restaurant, food processing, beauty salon, electronics

Product Details:

- 1, Material: Non woven PP
- 2, Size: 17.5x9.3cm, 17.5x9.5cm, as customers' requirements.
- 3, Color: Blue(in Stock), White, Green, pink, yellow, grey is are also available.
- 4, Style: Elastic earloop or with ties on.
- 5, Layer: Usually 3ply, 2ply and 4ply are also available for different requirements.

- 6, Nose bar: Single aluminum strip or double aluminium strip
- 7, Economic, light, Non-irritating, Comfortable
- 8, Prevent contamination from bacteria, dust, or coughing

Model Number: CHM1
Function: Surgical Face Mask
Product name: Disposable Face Mask
Material: 3 ply Non Woven fabric made in Germany.
Type: Antibacterial
Application: Surgical/Protection
Packing: 50pcs/Bag
Discription: Disposable Surgical Face Mask

Packaging & Delivery
Selling Units: Box of 50 masks.
Single package size: 17.5X9.5X1 cm

Masks are tested and approved to EN14683 standard.

Prices may change at any time.

Product and packaging pictures may be different from the delivered product.

Size of mask(Unexpanded) : 9.5cm x 17cm



