



BIOIDE-95 COMFORT PROTECTION

**€1,000.00 Original price was:
€1,000.00. ~~€49.00~~ Current price is:
€49.00.**



**STOCK CLEARANCE £0.05. ONLY £49 for
box with 1000pcs!**

Excellent 5 layer facemasks. Giveaway at
only £0.05

Box of 1000pcs £49.

Each individual packed in a protective
sleeve, 10 sleeves in a bag, 100 bags in the
box.

Immediate shipment from London
warehouse.

Bioide Pharmaceuticals Wholesale: Type-95
multi layer filtration. Comfortable to wear
and fits tightly around nose and mouth for
maximum protection.

Tested by SGS and TUV Nord. CE and FDA
certified. Packed in a protective sleeve.
Minimum Order: 1 Carton of 1000 pcs from
Bioide Warehouses in UK.

SKU: DCC95WS

Category: [Offers](#)

PRODUCT DESCRIPTION

Click for Detailed Description

CE and FDA certified - Quality Checks by SGS and TUV Nord.

FFP2 KN95 Disposable

Material: 5ply(non-woven+Hot air cotton+meltblown cloth+non-woven)

Color: White

Standard: China GB2626-2006 KN95 = N95 = FFP2

Application: Personal protection, PM2.5 protection

Filtering effect:95%

Use for hospital,laboratory,food industry,dust free workshop,electronic manufacturers.

Product Features

- 1.Multi layer filter structure protection
- 2.Ear belt / head-mounted
- 3.Comfortable and breathable
- 4.Good protective effect
- 5.Comfortable soft elastic band
- 6.Soft sponge nose pad
- 7.Aluminum adjustable nose clip

General information - Contact us for more details if needed:

Quality checked by several laboratories. CE certificates verified
Supplied from multiple factories. Photos can vary by manufacturer.

Packaging: Op-Bags 10 pcs / Box 50pcs / Carton 1000 pcs.

QTY in carton: 1000 pcs.

[Bioide-Type95_TechFile220520](#)

