



## **COVERALL - CORONA VIRUS DISPOSABLE COVERALL PROTECTION CLOTHING**

**€14.90**

High quality Corona Virus Disposable  
Coverall

Surgical hospital Protection Clothing

Size: m,l,xl

Weight: 20+20+25 gsm or

23+25/BFE99%+34 gsm

Usage: Indoor & Outdoor

Application: Isolation gown

Style: Antibacterial

Material: pp nonwoven+meltblown

nonwoven+pp nonwoven

Packing: 10pcs/box

---

**SKU:** N95Overall

**Categories:** [Bioide](#), [PersonProtect](#)

## **PRODUCT DESCRIPTION**

Corona Virus Disposable Coverall Surgical hospital Protection Clothing

Material SBPP/PE/PP+PE/SMS/SMMS/SF Microporous Material, 30~50g, or as your required Design

\*Breathable and waterproof

\*Acid-proof, Anti-Shrink, Anti-Static, Breathable, Dust-Proof, Greaseproof, Nontoxic, Plus Size, Quick Dry and so on.

\*Hood, waist, cuffs and ankles with elastic adjustment, A zip front and elastic cuffs and ankles

\*Tough and fiber free

\*Special for cleaning,restoration,industry,painting,laboratories,etc

\*It also can printed your specific Logo, different specifications are available!

\*High quality,wholesale price,fast delivery

Size XS, S, M, L, XL, XXL,XXXL

Color as picture or can be customized

Feature anti-shrink , Breathable, anti-wrinkle, eco-friendly, washable

